

# Insul-Eze®

---

## Frequently Asked Questions/ Troubleshooting

**Q: Can you suggest an easier method to use the Insul-Eze® with one hand?**

**A:** Gluing a magnet to the back of the Insul-Eze and applying it to your refrigerator or other magnetized surface will help you to hold the device steady while you load your insulin syringe. Velcro\* used in a similar manner may also be helpful.

**Q: Does the Insul-Eze contain latex?**

**A:** The Insul-Eze does not contain latex. For a list of materials please contact AMBIMEDINC.



**AMBIMEDINC**  
World Class Healthcare Products

908 Sir Francis Ave Capitola, CA 95010 • 831-475-1765 phone/fax  
info@ambimedinc.com • www.ambimedinc.com

Insul-Eze is a registered trademark of AMBIMEDINC. All rights reserved.  
\*Other trademarks are property of their respective holders.