


AMBIMEDINC Products



Inject-Ease

Injections Made Easy

Conquer the technique and anxiety of injections with the Inject-Ease, the device that makes injections easy. Simply place a loaded syringe in the Inject-Ease, place the tip against the skin, and press the button to automatically deliver the needle through the skin. You control the rate at which the medicine is injected.

- Usable with subcutaneous injection treatments, including diabetes and other conditions
- Special tip & technique- helps reduce the pain by masking the needle puncture
- One press of the button inserts the needle through the skin automatically and easily
- Part #IE400




Insul-eze™

Syringe Magnifier

One-piece syringe magnifier that holds syringe and insulin bottles in place for those with visual or manual dexterity problems.

- Aids in syringe loading for visually impaired
- Magnifies all syringe calibrations two times for easy reading
- Holds syringe and bottle firmly while drawing dosage
- Part #IZ6000



Insul-Cap™

Stores, Identifies Insulin Bottles, and Helps Load Syringes

Insulin bottle caps designed to aid in storing, identifying and making drawing insulin into a syringe easier.

- Allows one-handed syringe loading
- Fits most standard syringes and insulin bottles
- Keeps bottle tops sanitary
- Part #IC200S



Insul-Tray™

Insulin Bottle Tray

Small and sturdy Insulin bottle holder that keeps bottles organized and standing upright.

- Fits most standard insulin bottles
- Easy to put on, store and remove
- Part #ITR500S



Insul-Tote™

Insulated Carry Totes

Compact and durable all-weather insulated totes for transporting temperature sensitive medication, supplies, and meters.

- All-weather insulated zippered carry tote
- Transports temperature-sensitive medication
- Made of washable nylon with foam insulation



auto-Lancet

Lancing Device with Adjustable Tip

Re-usable quiet and comfortable lancing device that assists in obtaining capillary blood samples for testing.

- Five level adjustable tip allows variation of lancet penetration
- Pen shape available in standard (5") or mini (3.5") size
- Works with industry standard lancets
- Five year warranty



EZ-Lets™

Lancets

High-quality lancets compatible with industry standard lancing devices.

- Precision-ground tri-bevel point for clean, precise punctures
- Full range of sizes available — 21 gauge, 26 gauge, and 28 gauge



EZ-Lance®

The Ultimate Disposable Safety Lancet

The EZ-Lance® maximizes comfort, safety, and ease of use. Every feature has been designed with users in mind, from the extra-large tab for easy removal to the ergonomic thumb-groove for the utmost in comfort.

Safety

- Single-use, pre-cocked and not re-cockable
- Needle is covered at all times before and after lance
- Easy to hold; large tab is easy to remove



EZ-Vac™

Vacuum Alternate Site Testing Device

The EZ-Vac works with any lancing device to allow blood to be drawn from areas other than the fingertip, such as the forearm.

- Compact design — smallest, most portable vacuum Alternate Site Testing device
- Multi-location tip — Special tip works on alternate locations