

EZ-Lance[®]

The Ultimate Safety Lancet

The revolutionary new design of the **EZ-Lance[®]** maximizes comfort, safety, and ease of use. Every feature has been designed with users in mind, from the extra-large tab for easy removal to the ergonomic thumb-groove for the utmost in comfort. The high-speed, linear-tracking lance is highly advanced technology designed to help minimize pain. The internal breakaway bar is designed to prevent reuse, while the safety indicator window on the bottom ensures that you know whether the **EZ-Lance** is used or unused.

Don't settle for less than the best - demand the **EZ-Lance!**

From the people at **AMBIMEDINC**, makers of world-class medical device products.

Ideal for use in
hospitals and labs



Key Features

Safety

- Single-use to eliminate cross-contamination
- Not re-cockable - maximum safety feature prevents reuse
- Safety indicator shows clearly if device has been used
- No needle exposure for maximum safety; needle is covered at all times before and after lance

Comfort

- Ergonomic design - easy to hold; large tab is easy to remove
- High-speed lance - helps minimize pain
- Linear tracking - minimizes needle vibration for a cleaner lance
- Special bounce-back restraint - single lance only, no chance for a second lance

Ease of Use

- Advanced design - minimal user instructions
- Pre-cocked - simplifies use; no extra step needed

Key Benefits

- Four different versions available
- Maximum safety - impossible to re-cock, lance is always concealed except at the moment of lancing
- Ultimate ease of use - remove tab, apply to skin, and press button to fire
- Ergonomic design for maximum comfort
- Advanced lancing technology to minimize pain

Ordering Information

AMBIMEDINC Order Hotline:
1-831-475-1765
E-mail: info@ambimedinc.com

Please note that not all products are available in all markets. For full details, please contact your local reseller or AMBIMEDINC.

Special Orders

For special configurations and high volume orders, contact AMBIMEDINC or send e-mail to info@ambimedinc.com.

Specifications & Packaging Configurations

Yellow Model - Lite Flow

Color: White (body) and yellow (tab)
Needle size: 26 gauge
Point configuration: Tri-bevel
Depth of Penetration: 1.8 mm



PN/REF: ELY100 (100 ct. box)
UPC Code: 0-25294-04814-9

Green Model - Regular Flow

Color: White (body) and green (tab)
Needle size: 21 gauge
Point configuration: Tri-bevel
Depth of Penetration: 1.8 mm



PN/REF: ELG100 (100 ct. box)
UPC Code: 0-25294-04808-8

Orange Model - Regular Plus Flow

Color: White (body) and orange (tab)
Needle size: 21 gauge
Point configuration: Tri-bevel
Depth of Penetration: 2.2mm



PN/REF: ELO100 (100 ct. box)
UPC Code: 0-25294-04810-1

Red Model - Ultra Flow

Color: White (body) and red (tab)
Needle size: 21 gauge
Point configuration: Microblade
Depth of Penetration: 2.2mm



PN/REF: ELR100 (100 ct. box)
UPC Code: 0-25294-04812-5

All Models

Length (w/tab): 1.87 in. (47.49 mm)
Largest width: 0.79 in. (20.06 mm)
Height: 0.54 in. (13.71 mm)
HCPC Code (all models): A4259

EZ-Lance is a registered trademark of AMBIMEDINC,
U.S. Patent # 7,282,058
Canada Design Registration #105226
Registered Community Design # 0001139020001
International patents pending.
All Rights Reserved.

AMBIMEDINC
World Class Healthcare Products

908 Sir Francis Avenue
Capitola, CA • 95010, USA
Phone/Fax 1-831-475-1765
e-mail: info@ambimedinc.com
web: www.ambimedinc.com

Specifications listed are preliminary. Final specifications subject to change without notice.
© 2010, AMBIMEDINC, Inc.