

Company History

Since 1990, AMBIMEDINC has been providing high quality and innovative products to the clinical lab market. In 2009, AMBIMEDINC acquired the product line from Palco Labs, a company founded in 1979, as a leading diabetes product manufacturer of reusable and disposable devices for collecting capillary blood samples as well as injection aids for the diabetes and Urology impotence markets.

Mission

AMBIMEDINC strives to help people lead healthy and productive lives through the design and manufacture of high quality, innovative medical device products.

Values

- **Quality** products distinguish AMBIMEDINC in a competitive market. We strive to deliver products of superior quality through continuous improvement and manufacturing excellence.
- **Customer commitment** is the focus of our business. We measure our success through the continuing satisfaction of our customers.
- **Innovation** is the foundation of our business. We apply creative design and innovative engineering to craft new medical device products for the markets of today and tomorrow.
- **Employee satisfaction** is key to all of our company values. We strive to create a workplace based on team work and fairness, which promotes employee growth and participation.
- **Corporate citizenship** signifies our responsibility to the community and the environment as a company and as individuals.

Product Overview

AMBIMEDINC manufactures world class medical device products for distribution in more than 40 countries. Known for its quality standards and innovative designs, AMBIMEDINC focuses on medical device products designed to assist with blood sampling and injections.

Lancing



auto-Lancet™ "Classic"

The original adjustable tip. This reusable lancing device (in standard and mini sizes) assists in obtaining capillary blood samples for testing. The auto-Lancet is constructed in the USA of durable plastics for precision lancing day after day, year after year. It comes with a 5 year warranty.

auto-Lancet™ "Custom"

Modern version of the original adjustable tip. AMBIMEDINC's newest auto-Lancet delivers fully customizable colors to match any meter. The auto-Lancet "Custom" comes in a new, compact size and offers technology for faster lancing and improved linear-tracking. Comes with a 5 year warranty.

Company Fact Sheet

Lancing



EZ-Lets™

Experience the least amount of pain with AMBIMEDINC's new blood lancets. EZ-Lets are premium-quality, precision-engineered sterile blood lancets. The finely-ground needle edges reduce discomfort and encourage the healing process. Offers easy-to-use, twist-off cap. EZ-Lets fit most reusable finger lancing devices, and are available in three sizes: 21 gauge, 26 gauge and 30 gauge.

Assistive



Inject-Ease®

Injections made easy. Conquer the anxiety of injections with the Inject-Ease, the device that makes injections easy. Simply place your loaded syringe in the Inject-Ease, place the tip against your skin, and press the button to automatically deliver the needle through the skin. You control the rate of injection. The special tip is designed to help reduce the pain of injection. Perfect for many subcutaneous injections including those for diabetes, arthritis, psoriasis, asthma, infertility, vaccinations and others.

Lancing



EZ-Lance™

The Ultimate Disposable Safety Lancet. The design of the EZ-Lance® maximizes comfort, safety, and ease of use. Every feature has been designed with users in mind, from the extra-large tab for easy removal to the ergonomic thumb-groove for the utmost in comfort. The high-speed, in-line needle strike minimizes the perception of pain. Safety features eliminate the possibility of reuse, and show whether the EZ-Lance is used. Four lancets to choose from.

Patents

AMBIMEDINC has a number of patented innovations in its medical devices, and continually innovates to create new patents.

Trademarks

EZ-Lets, Inject-Ease and EZ-Lance are registered trademarks of AMBIMEDINC. All rights reserved.



908 Sir Francis Ave Capitola, CA 95010 • 831-475-1765 phone/fax
info@ambimedinc.com • www.ambimedinc.com