

# auto-Lancet™

## Enduring Quality

### The Original Adjustable Tip

The incomparable **auto-Lancet™** is the lancing device that has set the standard for quality, comfort, and adjustability since 1983. This re-usable lancing device assists in obtaining capillary blood samples for testing. The **auto-Lancet** is constructed in the USA of durable plastics for precision lancing day after day, year after year.

The **auto-Lancet** - quality you can count on.  
From the people at AMBIMEDINC,  
makers of world-class medical device  
products and leaders in  
lancing technology.



### Key Features/Benefits

- Adjustable tip - Allows precise customization of lancet penetration
- High quality construction - Five (5) year warranty
- Made of durable plastics for years of use
- Convenient pen shape available in standard or mini size
- Made in the USA

### Key Selling Points

- Available in standard or mini size
- Works with industry standard lancets
- #1 lancing device chosen by glucose meter manufacturers worldwide

## Ordering Information

AMBIMEDINC Order Hotline:

1-831-475-1765

E-mail: [info@ambimedinc.com](mailto:info@ambimedinc.com)

Please note that not all products are available in all markets. For full details, please contact your local reseller or AMBIMEDINC.

## Special Orders

For special configurations and high volume orders, contact AMBIMEDINC or send e-mail to [info@ambimedinc.com](mailto:info@ambimedinc.com).

auto-Lancet is a trademark of AMBIMEDINC.  
All rights reserved.

\*Other trademarks are property of their respective holders.

**AMBIMEDINC**  
World Class Healthcare Products

908 Sir Francis Avenue  
Capitola, CA • 95010, USA  
Phone/Fax 1-831-475-1765  
e-mail: [info@ambimedinc.com](mailto:info@ambimedinc.com)  
web: [www.ambimedinc.com](http://www.ambimedinc.com)

Specifications are subject to change without notice.  
© 2010 AMBIMEDINC. All rights reserved.

## Specifications



auto-Lancet Standard (actual size)

### Standard auto-Lancet Specifications:

Length: 4.905 in. (12.45 cm)

Barrel Width: 0.475 in. (1.206 cm)



auto-Lancet Mini (actual size)

### Mini auto-Lancet Specifications:

Length: 3.595 in. (9.13 cm)

Barrel Width: 0.470 in. (1.19 cm)

### Lancet Compatibility:

Industry standard, including AMBIMEDINC EZ-Lets, Ultra-Fine\* II, Monolet\*, Unilet\* GP, Microlet\*, OneTouch\*, Ultra TLC\*, SoftTouch\* [Roche Australia], SoftTouch\* [Roche Accu-Chek], Gentle-Let\*, CleanLet\*, Omni-can\*, Qualitilet\*, and Freestyle\*.

**HCPC Code:** A4259

### Warranty:

AMBIMEDINC warrants that each new auto-Lancet shall be free of all manufacturing defects and guarantees normal performance for five (5) years.