

FOR IMMEDIATE RELEASE

AMBIMEDINC MODIFIES INJECT-EASE TO MAKE IT COMPATIBLE WITH TERUMO & MONOJECT INSULIN SYRINGES FOR THE DIABETIC MARKET

Capitola, California – June 2, 2014 – AmbiMedInc, Manufacturer of Inject-Ease, a device that makes injections easy and comfortable, announced modifications to its product making it not only compatible with BD insulin syringes as originally designed, but also compatible with many other insulin syringes such as Terumo, Monoject and others that have similar dimensions.

The Inject-Ease is an auto injector, usable for subcutaneous self-injections and is designed to make injections easy and comfortable. The syringe is simply loaded in the Inject-Ease then the tip of the device placed against the skin and rotated to create less feeling when the needle goes into the skin. A button is then pressed to automatically deliver the needle through the skin and the user inserts the medicine.

When using the Inject-Ease, one does not see the needle thus reducing anxiety. Pain is minimized by the twisting of the tip on the skin and when the button is pushed, the needle is inserted into the skin rather than the patient doing it.

The Inject-Ease is compatible with insulin syringes from BD, Terumo, Monoject and many others syringes that have similar dimensions. It is ideal for anyone who has a fear of needles and designed for one hand use allowing access to hard to reach injection sites. It is also a supreme learning tool for the new patient.

The Inject-Ease comes with 3 adapter tips (30 unit for 0.3 ml, 1/3cc syringe; 50 unit for .5 ml, 1/2cc syringe; 100 unit for 1.0 ml, 1cc syringe), 2 spacer rings (1/8 in., 0.317 cm; 1/4 in., 0.635 cm) and a clear plunger cap for preloading medications.

About AmbiMedInc

AmbiMedInc is a world class medical device manufacturer of reusable and disposable devices for collecting capillary blood samples as well as injection aids. Its products are distributed in more than 80 countries and it is known for its high quality standards and innovative designs.

To get more information about the Inject-Ease, please contact AmbiMedInc at info@ambimedinc.com or 831-475-1765 or access www.ambimedinc.com
http://www.ambimedinc.com/section_products/injectease.html

Visit AmbiMedInc at the American Association of Diabetes Educators annual meeting in Orlando, August 6-9, 2014, **Booth#856**.