


Insul-Cap



Frequently Asked Questions/ Troubleshooting

Q: What size syringe works with the Insul-Cap?

A: The Insul-Cap works with all size BD insulin syringes.

Q: What size vial works with with the Insul-Cap?

A: The Insul-Cap fits Eli Lilly*, Novo Nordisk*, and Aventis Lantus* insulin bottles.

Q: Do the short BD* syringes work with the Insul-Cap?

A: Yes, the Insul-Cap works with the short syringes.

Q: Does the Insul-Cap contain latex?

A: The Insul-Cap does not contain latex. For a list of materials please contact AMBIMEDINC.

Q: How do I clean the Insul-Cap?

A: Submerge the Insul-Cap in a soap and water solution, rinse and air dry.

