



908 Sir Francis Avenue, Capitola, CA 95010 USA  
Phone/Fax 831 475-1765  
[www.ambiMEDinc.com](http://www.ambiMEDinc.com)

Media Contact:  
Ellen Brady

FOR IMMEDIATE RELEASE

## **AMBIMEDINC RELEASES INSUL-EZE - A NEW AND IMPROVED SYRINGE MAGNIFIER ENABLING EASY SYRINGE LOADING FOR THOSE WITH POOR VISION AND DEXTERITY**

**Capitola, California** – October 7, 2013 – AmbiMedInc, Designer and Manufacturer of World Class Medical Devices released the new and improved **Insul-eze** syringe magnifier at the American Association of Diabetes Educators Annual Meeting in Philadelphia.

The **Insul-eze** received astounding praise from the Educator community for its sturdy one-piece design that holds the syringe and insulin bottle in place for those with visual and/or manual dexterity problems while magnifying the calibrations over two times for easy reading.

Recent modifications allow all size insulin bottles to securely snap into place while a channel on the underside allow all sizes of insulin syringes to easily fit and effortlessly align the needle into the bottle top. A new sturdy base prevents rolling, creates stability and allows the magnifier, syringe and bottle to be secure for easy loading, viewing and ease of use. The product can also be affixed to a counter or refrigerator with Velcro or magnet to allow one-handed loading.

### **About AmbiMedInc**

AmbiMedInc is a U.S. manufacturer of world class medical device products for distribution in more than 40 countries and is known for its quality standards and innovative medical products. As a leading diabetes product manufacturer of reusable and disposable devices for collecting capillary blood samples as well as injection aids for the diabetes, AmbiMedInc continues to strive to bring new and helpful products to the diabetic community.

To get more information about the Insul-eze, please contact AmbiMedInc at [info@ambimedinc.com](mailto:info@ambimedinc.com) or 831-475-1765 or access [www.ambimedinc.com](http://www.ambimedinc.com).