

auto-Lancet™

Enduring Quality

The Original Adjustable Tip

The incomparable **auto-Lancet™** is the lancing device that has set the standard for quality, comfort, and adjustability since 1983. This re-usable lancing device assists in obtaining capillary blood samples for testing. The **auto-Lancet** is constructed in the USA of durable plastics for precision lancing day after day, year after year.

The **auto-Lancet** - quality you can count on.
From the people at AMBIMEDINC,
makers of world-class medical device
products and leaders in
lancing technology.



Key Features/Benefits

- Adjustable tip - Allows precise customization of lancet penetration
- High quality construction - Five (5) year warranty
- Made of durable plastics for years of use
- Convenient pen shape available in standard or mini size
- Made in the USA

Key Selling Points

- Available in standard or mini size
- Works with industry standard lancets
- #1 lancing device chosen by glucose meter manufacturers worldwide

Ordering Information

AMBIMEDINC Order Hotline:
1-831-475-1765
E-mail: info@ambimedinc.com

Please note that not all products are available in all markets. For full details, please contact your local reseller or AMBIMEDINC.

Special Orders

For special configurations and high volume orders, contact AMBIMEDINC or send e-mail to info@ambimedinc.com.

auto-Lancet is a trademark of AMBIMEDINC.
All rights reserved.

*Other trademarks are property of their respective holders.

AMBIMEDINC
World Class Healthcare Products

908 Sir Francis Avenue
Capitola, CA • 95010, USA
Phone/Fax 1-831-475-1765
e-mail: info@ambimedinc.com
web: www.ambimedinc.com

Specifications are subject to change without notice.
© 2010 AMBIMEDINC. All rights reserved.

Specifications



auto-Lancet Standard (actual size)

Standard auto-Lancet Specifications:

Length: 4.905 in. (12.45 cm)
Barrel Width: 0.475 in. (1.206 cm)



auto-Lancet Mini (actual size)

Mini auto-Lancet Specifications:

Length: 3.595 in. (9.13 cm)
Barrel Width: 0.470 in. (1.19 cm)

Lancet Compatibility:

Industry standard, including AMBIMEDINC EZ-Lets, Ultra-Fine* II, Monolet*, Unilet* GP, Microlet*, OneTouch*, Ultra TLC*, SoftTouch* [Roche Australia], SoftTouch* [Roche Accu-Chek], Gentle-Let*, CleanLet*, Omni-can*, Qualitilet*, and Freestyle*.

HCPC Code: A4259

Warranty:

AMBIMEDINC warrants that each new auto-Lancet shall be free of all manufacturing defects and guarantees normal performance for five (5) years.