


AMBIMEDINC PRODUCTS



Inject-Ease® *Injections Made Easy*

Conquer the technique and anxiety of injections with the Inject-Ease, the device that makes injections easy. Simply place a loaded syringe in the Inject-Ease, place the tip against the skin, and press the button to automatically deliver the needle through the skin. You control the rate at which the medicine is injected.

- Usable with subcutaneous injection treatments, including diabetes and other conditions
- Special tip & technique- helps reduce the pain by masking the needle puncture
- One press of the button inserts the needle through the skin automatically and easily
- Part #IE400



Insul-eze® *Syringe Magnifier*

One-piece syringe magnifier that holds syringe and insulin bottles in place for those with visual or manual dexterity problems.

- Aids in syringe loading for visually impaired
- Magnifies all syringe calibrations two times for easy reading
- Holds syringe and bottle firmly while drawing dosage
- Part #IZ600



Insul-eze®SPM *Magnifier for Large Syringe and Pen Dial Injectors*

Clip-on magnifier for syringes larger than 1ml, and Pen Injectors. This easy-to-use magnifier allows those with visual impairments to easily see the level of medication on the syringe body and to see the dose dial on pen injectors.

- Magnifies syringe calibrations two times for easy reading
- Increases viewing area on the dose dial of pen injectors
- Fits the barrel of 3ml, 5ml and 10ml syringes, and most pen injectors
- Easy-Squeeze mechanism makes it simple to put on and take off
- Part #IZ600SPM



Insul-Cap® *Stores, Identifies Bottles, and Helps Load Syringes*

Insulin bottle caps designed to aid in storing, identifying and making drawing insulin into a syringe easier.

- Allows one-handed syringe loading
- Fits most standard syringes and insulin bottles
- Keeps bottle tops sanitary
- Part #IC200S



Hide-a-Needle® Needle-Hider for Pen Injectors

A simple device that locks onto Pen Injectors, hiding the needle when administering subcutaneous injections. The Hide-A-Needle eases the anxiety of injections, and its special design helps reduce the pain by masking the needle puncture.

- Works with most pen devices
- Easy to put on and off
- Hides the pen injector needle
- Part #HAN300



Ambi-Tray® Insulin Bottle Tray

Small and sturdy Insulin bottle holder that keeps bottles organized and standing upright.

- Fits most standard insulin bottles
- Easy to put on, store and remove
- Part #ITR500S



Insul-Tote® Insulated Carry Totes

Compact and durable all-weather insulated totes for transporting temperature sensitive medication, supplies, and meters.

- All-weather insulated totes
- Transports temperature-sensitive medication
- Made of washable nylon with foam insulation
- Includes handy pockets for storing supplies and reusable frozen gel pack
- 5 Models to choose from:
 - Pen-Tote - Part #IT500PTD, for shorter outings and the most discreet travel. Comes with 3oz gel pack. Size 4in x 7in x 1in (1.8cm x 10.2cm x 2.5cm)
 - Hip-Pack - Part #IT500JR, waist pack with adjustable strap for convenient hands-free carrying. Comes with 5oz gel pack. Size 6in x 4.5in x 3in (15.2cm x 11.4cm x 7.6cm)
 - Top-Opening - Part #IT500T, bag style with shoulder strap for easy carrying. Comes with 9 oz gel pack. Size 7in x 7in x 3in (17.8cm x 17.8cm x 7.6cm)
 - Back Opening - Part #IT500B, opens out flat for quick viewing of stored items. Comes with 9 oz gel pack. Size 7in x 7in x 3in (17.8cm x 17.8cm x 7.6cm)
 - Back Opening Ultimate - Part #IT500U, wider in size with extra pockets for more storage. Opens out flat for quick viewing of stored items. Comes with 9oz gel pack. Size 7in x 7in x 4in (17.8cm x 17.8cm x 10.2cm)



ambi-Lancet® Lancing Device with Adjustable Tip

Re-usable quiet and comfortable lancing device that assists in obtaining capillary blood samples for testing.

- Five level adjustable tip allows variation of lancet penetration
- Works with industry standard lancets
- 3 Models to choose from:
 - ambi-Lancet Standard - Part #ALI00K23AESP, length 4.5 in (11.43 cm)
 - ambi-Lancet Mini - Part #ALI00MK23AESP, length 3.75 in (9.52 cm)
 - ambi-Lancet Eject - Part #ALI00N23AESP, length 4 in (10.16 cm)

AMBIMEDINC
World Class Healthcare Products