

ambi-Lancet®

Enduring Quality

Lancing Devices with Adjustable Tip

Re-usable lancing device assists in obtaining capillary blood samples for testing. The **ambi-Lancet®** is constructed of durable plastics for precision lancing day after day, year after year.

The **ambi-Lancet®**—quality you can count on.

From the people at AMBIMEDINC, makers of world-class medical device products and leaders in lancing technology.



Key Features/Benefits

- Adjustable tip - Allows precise customization of lancet penetration
- High quality construction - One(1) year warranty
- Convenient pen shape

Key Selling Points

- Available in three configurations; standard, mini and eject
- Works with industry standard lancets

Specifications

Ordering Information

AMBIMEDINC Order Hotline:
1-831-475-1765
Email: info@ambimedinc.com

Please note that not all products are available in all markets. For full details, please contact your local reseller or AMBIMEDINC.

Special Orders

For special configurations or high volume orders, contact AMBIMEDINC or send an email to info@ambimedinc.com

ambi-Lancet is a trademark of AMBIMEDINC.
All rights reserved.

*Other trademarks are property of their respective holders.

AMBIMEDINC
World Class Healthcare Products

Capitola, CA • 95010, USA
Phone 1-831-475-1765
Email: info@ambimedinc.com
web: www.ambimedinc.com

Specifications are subject to change without notice.
© 2018 AMBIMEDINC. All rights reserved.



ambi-Lancet Standard (actual size)

ambi-Lancet Standard Specifications:

Length: 4.5 in. (11.43 cm)
Barrel width: 0.55 in. (1.34 cm)
Part #AL100K23AESP



ambi-Lancet Mini (actual size)

ambi-Lancet Mini Specifications:

Length: 3.75 in. (9.52 cm)
Barrel width: 0.55 in. (1.34 cm)
Part #AL100MK23AESP



ambi-Lancet Eject (actual size)

ambi-Lancet Eject Specifications:

Length: 4 in. (10.16 cm)
Barrel width: 0.63 in. (1.60 cm)
Part #AL100N23AESP

Lancet Compatibility:

Industry standard, including AMBIMEDINC EZ-Lets, Ultra-Fine* II, Monolet*, Unilet* GP, Microlet*, OneTouch*, Ultra TLC*, SoftTouch* [Roche Australia], SoftTouch* [Roche Accu-Chek], Gentle-Let*, CleanLet*, Omni-can*, Qualitilet*, and Freestyle*.

HCPC Code:A4259

Warranty:

AMBIMEDINC warrants that each new auto-Lancet shall be free of all manufacturing defects and guarantees normal performance for 1 year.