

Insul-eze SPM

Magnifier for Pen-Injector Dial and Large Syringes

Frequently Asked Questions/Troubleshooting



Q: What syringes work with the Insul-eze SPM?

A: The Insul-eze SPM fits the barrel of most 3ml, 5ml, 10ml syringes and 200-unit syringes.

Q: What Pen-Injectors work with the Insul-eze SPM?

A: The Insul-eze SPM fits on the dial of most Pen-Injectors.

Q: How much does the Insul-eze SPM magnify?

A: The Insul-eze SPM magnifies the syringe calibrations 2 times for easier reading.

Q: Does the Insul-eze SPM contain latex?

A: The Insul-eze SPM does not contain latex. For a list of materials please contact AmbiMedInc.